

Automated Pipette Tip Washer

Wash & Dry Throughput: 6-8 Tip Racks / Hour
Dimensions: 16W x 24H x 16D inches

Pipette Tip Compatibility:

- Any volume range 10 µl- 5mL
- Any tip rack format; 96, 384, etc.
- Handheld Manual Pipette Tips
- Automation Pipette Tips
- Clear Pipette Tips
- Black Conductive Tips



Liquid Handling Automation Integratable

Rethinking Laboratory Waste



TIPNOVUS™
mini



Wash features: ▾

- ▣ **Four cleaning reagent feeds:** Enables the users to wash and clean consumables with the cleaning reagent of their choice
- ▣ **High Pressurized wash:** Removes any large contaminants from pipette tip surfaces
- ▣ **Ultrasonic Cleaning:** Removes small contaminants and particles from pipette tip surfaces
- ▣ **Tip Agitation:** Improves the wash quality and overall cleaning process
- ▣ **UV sanitation:** Kills 99.9% of infectious diseases from tip surfaces and makes tips DNA, RNA free

LEAF ▾

An intuitive, user friendly graphical software, LEAF™, allows users to program their own custom wash protocol for the TipNovus. LEAF™ software will enable the users to select the appropriate cleaning solutions, number of wash cycles, sonication and UV sanitation features for the custom wash protocol.



Dry features: ▾

- ▣ **Tip Agitation:** Remove large droplets and improve overall drying process
- ▣ **Heat-Fan Dry:** Automated temperature Control Heat & Fan for optimal drying
- ▣ **Optional UV Sanitation**



GRENOVA

Phone: 804-325-1354 - Fax: 804-299-5208 - sales@grenovasolutions.com - www.grenovasolutions.com